

**NEW**



# **HANDLER**

**PRO TECH<sup>®</sup>**  
**ELITE**

## **55VF Vacuum Former**



**The 55VF is Handler's latest addition to a robust portfolio of dental lab products.**

This is an economical way to process:

- Bleaching trays
- Mouthguards
- Implant splints
- Temporary splints
- Custom impression trays
- Denture bases
- Copings
- Temporary crowns and bridges
- Orthodontic appliances

The 55VF can accommodate nearly all dental vacuum forming applications.

- Machine weight: 25.5 lbs.
- Machine size: 10-11/16"W x 17-13/16"H x 10-5/8"D

### **Features:**

- Sturdy stainless steel, powder coating and aluminum construction.
- Reversible vacuum tray to manage multiple applications.
- Container of stainless steel beads.
- Easy to operate design.
- 750 watt symmetrical heating element heats the material quickly and evenly.
- Heat insulation design is safe for the user, and cooler to the touch on the top cover.
- Ergonomic handle lift.
- The no twist top makes the operation smoother and easier.
- When you bring down the heated material, the vacuum activates automatically.
- Stainless steel bearing blocks on the head assembly allow for smooth up and down motion.
- Spring Action handle for easy open and close of the frame.
- Uses all types of 5" x 5" resin sheets.